

## Polysome fractionation with RNA purification

Polysome fractionation with RNA purification  
Johanna Schott, Georg Stoecklin lab, 2013

### Material:

Sucrose

Beckman Centrifuge Tubes, 11 x 60 mm (328874)

Cyclohexemide (CHX)

1 x PBS (Phosphate buffered saline)

Ultracentrifuge

SW60 rotor

Fractionator (Teledyne ISCO UA-6)

Phenol : Chloroform : Isamylalcohol (PCI) 25 : 24 : 1 (AppliChem, A0944)

Isopropanol

GlycoBlue (Ambion, AM9515)

Ethanol

RNAse-free water

### **2 x Polysome buffer**

30mM Tris-HCl pH7.4

30mM MgCl<sub>2</sub>

600mM NaCl

### **Polysome lysis buffer (for 10 ml)**

2 x Polysome buffer

5 ml

Water

4.8 ml

Triton-X-100

100  $\mu$ l

CHX (10 mg/ml)

100  $\mu$ l

complete protease inhibitor  
(Roche cOmplete Mini,  
EDTA-free)

1 tablet

RNasin (Promega)

50  $\mu$ l

$\beta$ -mercaptoethanol

10  $\mu$ l

### **Solution II (for 100 ml)**

1 M Tris-HCl pH 7.5

1 ml

5 M NaCl

7 ml

0.5 M EDTA

2 ml

10% SDS

10 ml

Urea

42 g

water

ad 100 ml

### Preparation of gradients:

> dissolve sucrose in 1 x Polysome buffer: 17.5%, 25.6%, 33.8%, 41.9% and 50%

> start with 790  $\mu$ l of 50% sucrose per tube

> freeze at -80°C

> continue with 41.9% sucrose, freeze etc.

> store gradients at -80°C

### Protocol:

Georg Stoecklin  
Johanna Schott

A200

g.stoecklin@dkfz.de

## Polysome fractionation with RNA purification

### 1. On the day before:

- > seed cells: they should not be confluent, otherwise they might reduce translation;  
for adherent cells, a 10 cm dish per gradient is a good starting point.
- > thaw gradients in the cold room overnight (4°C)
- > prepare 60 % sucrose

### 2. Polysome fractionation

- > add 100 µg/ml cycloheximide (CHX) to all dishes
  - > incubate for 5 min at RT
  - > optional: wash once with ice-cold 1 x PBS containing 100 µg/ml CHX
  - > depending on the cell type:  
scrape cells off in ice-cold 1 x PBS + 100 µg/ml CHX and pellet by centrifugation at 4°C  
or: lyse directly in the dish
  - > lyse in 350 µl ice-cold polysome lysis buffer per sample
  - > tumble for 10 min in the cold room
  - > centrifuge at 9'300 g, 4°C for 10 min (to remove the nuclei)
  - > take the supernatant into a fresh eppendorf tube without touching the nuclear pellet
  - > take 30 µl for isolation of total cytoplasmic RNA (into 270 µl polysome lysis buffer + 300 µl Solution II )
  - > load about 300 µl lysate per gradient (fill up with lysis buffer if necessary; work in the cold room)
  - > centrifuge for 2.5 hrs, 35'000 rpm at 4°C in a SW60 rotor (Acceleration: 7, Deceleration: 1)
  - > collect 16 fractions of about 250 - 300 µl into eppendorf tubes that already contain 300 µl of Sol II.
- We elute the gradients using an ISCO UA-6 gradient fractionator, which monitors UV absorbtion at 254 nm. Typical settings for fractionation are 14 sec per fraction at 50% pump speed.

### 3. RNA purification

- > add 300 µl Phenol-Chloroform-Isoamylalcohol (PCI, 25:24:1) to each fraction (under the fume hood!) and vortex
- > heat the samples for 10 min at 65°C (heating block)
- > open and close the eppendorf tubes (under the fume hood) to release the pressure
- > centrifuge for 20 min at 17'000 g at RT
- > transfer the aqueous (upper) phase (about 500 µl) to a fresh tube and add 1 µl GlycoBlue to each sample. GlycoBlue contains glycogen (to assist RNA precipitation) and a blue dye (so that you can see the RNA pellet after precipitation).

## Polysome fractionation with RNA purification

- > add 600  $\mu$ l Isopropanol
- > vortex
- > precipitate over night at -20°C
- > vortex briefly
- > centrifuge for 20 min at full speed and 4°C
- > take off and discard the supernatant (carefully with a pipet, do not disturb the pellet)
- > wash once with 800  $\mu$ l cold 70% ethanol
- > centrifuge for 10 min at full speed and 4°C
- > take off the supernatant again with equal care
- > add 20  $\mu$ l RNase-free water
- > dissolve the pellet by heating at 65°C for 10 min and vortexing

## 4. mRNA quantification

- > measure mRNA of interest in the different fractions using qPCR or Northern blot analysis. You can also pool fractions to reduce the number of samples.